
## **CytoSorbents**<sub>TM</sub>

# **CytoSorb Therapy**

# Decision support for ECMO/ECLS patients



**CytoSorb Therapy - REGAIN CONTROL** 

### **CytoSorb Therapy - Initiation**

## **CytoSorb Therapy - Continuation**



considering CytoSorb as an adjunctive therapy

Adequate causative therapy and technical conditions of ECMO therapy (catheter position, flow rate etc.) are the base for therapeutic success

#### Sufficient stabilization

#### Stop CytoSorb therapy

#### Insufficient stabilization

Only minor reduction of

> Consider new adsorber after 12 hrs. (early change) or at the latest after 24 hrs.



Re-evaluate every 12 to 24 hrs.

despite 2-3 adsorbers

#### CytoSorb therapy appears futile

# Potential Indications for CytoSorb Therapy in ECMO/ECLS Patients:

The following clinical conditions are characterized or aggravated by hyperinflammation, often with deterioration and shock. Cytosorb therapy may therefore be considered in addition to standard of care and treatment of the underlying cause in:

- Cardiogenic shock
- ECPR
- Bridge to VAD surgery
- ARDS with high vasopressor demand
- Post cardiotomy syndrome
- Infective endocarditis
- Septic shock
- Liver failure (removal of bilirubin)
- Rhabdomyolysis (removal of myoglobin)

#### CytoSorbents Europe GmbH

#### CytoSorbents Switzerland GmbH

Müggelseedamm 131 c 12587 Berlin | Germany V

T +49 30 65 49 91 45 F +49 30 65 49 91 46 support@cytosorbents.com c/o MGM GmbH Wielandstrasse 5 | 4153 Reinach BL | Switzerland

T +41 61 71 37 37 8 F +41 61 71 37 37 9 support@cytosorbents.com This decision guidance is non-binding and cannot replace the therapy decisions of the treating physician, who is in all cases responsible for the development and implementation of an adequate diagnostic and therapeutic plan for each individual patient.

The clinical and preclinical data and results obtained with the CytoSorb adsorber are not transferable to other products. CytoSorb should only be administered by personnel who have been properly trained in administration of extracorporeal therapies. CytoSorb is not available for commercial sale in USA.

CytoSorb and CytoSorbents are trademarks of the CytoSorbents Corporation, USA. © Copyright 2019, CytoSorbents Europe GmbH. All rights reserved. B1088R02EN2019